




Short Introduction to Static Light Scattering



**Static Light Scattering using the
ALV/CGS-3 Compact Goniometer System**

ALV-GmbH Laser Vertriebsgesellschaft mbH



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

Content

- **1.0. Static Light Scattering (SLS)**
- **1.1. Brief theory of Light Scattering**
- **2.0. Performing a SLS measurement**
- **2.1. Performing a SLS measurement with the ALV-Correlator Software**
- **3.0. Performing a SLS measurement with the ALV-Fit & Plot Software**




1.0. Static Light Scattering (SLS)

- a classical method for the determination of **Molecular Mass**, shape and (possible) interaction (**Radius of Gyration**, **Second Virialcoefficient**) as well as superstructure (**Static Structure Factor**) of particles in solution
- is used (among others) for the determination of the **Aggregation Number** (number of molecules present in a structure) of micelles or polymer coils
- requires an **angular dependent** measurement of the **scattered light**
- SLS measurements are (in conclusion) limited to the analysis of **monomodal systems**, as the dependency of the scattering intensity increases with the sixth power of the particle diameter, so larger contributions (like e.g. aggregates) would highly interfere the mean intensity result from the particles
- therefore very clean and almost **dust free measurement conditions** are required



1.1. Brief theory of Light Scattering

Light Scattering is caused by the **interaction of electrons of an isotropic particle/molecule with an influencing alternating electromagnetic field** (light wave), which sets the electrons in oscillation. Hereby, dipoles with orientation towards the polarization plane of the incident light are induced (Hertz dipoles), which consequently emit a corresponding light field of the same frequency, that shows on time average a characteristic radiation pattern.





Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The distance and angular dependency of such a Hertz dipole is given as:

$$\frac{I_S(r, \theta)}{I_0} = \frac{16 \pi^4 \alpha^2 \sin^2 \psi}{\lambda^4 r^2}$$

I_S = Scattering intensity

I_0 = Intensity of the incident light


Θ = Scattering angle

α = Polarizability

ψ = Angle between dipole axis and direction of observer (detector)

r = distance to observer (detector)

λ = wavelength of the incident light



The scattering vector q derives from:

$$|q| = \frac{4\pi n_d \sin(\theta / 2)}{\lambda}$$



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The **fluctuation theory of Einstein** solves the contradiction that in a homogeneous system, scattered light can be detected in a homogeneous system despite the fact that in a large scattered volume there always exist corresponding partial volumes which produce a 180° phase-shifted scattering and thus destroy all scattered light destructively. The reason for this is the spatial and temporal fluctuation of the polarizability. This fluctuation in a partial volume dV can be attributed to a fluctuation of the concentration.

As an important factor here the refractive index increment $\frac{dn}{dC}$

i.e. the change in the refractive index (dn) with the concentration (dC) at a constant temperature, is taken into account.

For linear polarized light and a particle size of $\leq \lambda/20$ this results in

$$I_s = I_0 \frac{4\pi^2 \left(\frac{dn}{dC} \right)^2 M^2}{\lambda^4 r^2 A^2}$$

A = Avogadro constant

M = Molecular weight of the particles

taking into account the refractive index of the solvent n_0 (for correcting the vacuum-related wavelength of the incident light) and the particle number N per volume

$$N = \frac{CA}{M}$$

A = Avogadro constant

M = Molecular weight of the particles

C = Particle concentration



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

One obtains for the scattering intensity per (detection) volume unit

$$I_S = I_0 \frac{4\pi^2 n_0^2 \left(\frac{dn}{dC}\right)^2 MC}{A\lambda^4 r^2}$$

with the so called **Rayleigh Ratio** R_g

$$R_g = \frac{I_S}{I_0} \cdot \frac{r^2}{2}$$

In which the scattering intensity must take account of the fact that the detector, as a function of the scattering angle θ , in each case considers a different scatter volume (for this purpose, the measured scatter intensity $I_{Messung}$ according to

$$I_S = I_{Messung} \cdot \sin \vartheta$$

Correspondingly converted to the surface-normalized scattering intensity), and the introduction of an optical constant K

$$K = \frac{2\pi^2 n_0^2 \left(\frac{dn}{dC}\right)^2}{A\lambda^4}$$



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

the **fundamental equation of light scattering** is obtained

$$\frac{KC}{R_9} = \frac{1}{M} + 2A_2c +$$


A_2 = second virial coefficient

It applies that A_2 is negative for attractive interaction and A_2 is positive for repulsive interactions

When applying Kc / R for different concentrations against the mass concentration c and subsequent extrapolation to small concentrations, the axis section (intercept) yields the inverse molar mass $1/M$. From the slope of the extrapolation line $q \rightarrow 0$, derived from the extrapolations of the concentrations to $q \rightarrow 0$ the value for the second virial coefficient A_2 is obtained



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System





Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

For particles larger than $\lambda / 20$ intramolecular interferences occur due to the presence of several scattering centers, i.e. light from different partial regions of the particle is in a fixed phase relation to one another. However, the interference is dependent on the measured scattering angle, and for two scattering centers it can be shown that the phase shift j is related with the scatter vector q

$$|q| = \frac{4\pi n_d \sin(\theta/2)}{\lambda}$$

and the distance r of the two centers

$$\varphi = -\vec{q}\vec{r}$$



As a result, more light is scattered into forward direction (0°) than into back direction (180°). Thus, e.g. the measured scattering intensity of a larger particle ($> \lambda / 20$) at $q = 45^\circ$ is higher than at $q = 135^\circ$, although the Rayleigh theory predicts the same intensity.



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System


The intramolecular interferences depend essentially on the phase shift so that the interference is influenced by the change in the scattering angle and thus q . Also, by changing the wavelength λ of the incident light, q can be varied. For larger wavelengths, q becomes smaller, i.e. in order to make smaller sub-regions of the to be measured particle (e.g., polymer or nanoparticles) visible, and as a result the distances r of the interfering scattering centers seem to be reduced, the scattering vector q must be increased with a constant phase shift.

Particle $> \lambda / 20$
multiple scattering centres





Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



q-range	Particle structure	Particle characteristics / properties
$qr \ll 1$	whole coil as particle	Molecular weight, Radius of Gyration
$qr < 1$	Particle shape	Arrangement of substructures
$qr \approx 1$	Details of the particle shape	Anisotropy, axis lengths
$qr > 1$	Coil structure	Chain structure (rods, helix structure, etc.)
$qr \gg 1$	Chain segments	Tacticity, conformation

The intramolecular interferences sum up the amplitudes of the electromagnetic field, which are induced by the individual scattering centers. The mass-dependent Rayleigh ratio for the infinitesimal small particles is multiplied by the form factor $P(\vartheta)$.

$$R_{\vartheta}(\text{reales Teilchen}) = P(\vartheta) \cdot R_{\vartheta}$$

In general, $P(\theta)$ is averaged over all possible spatial orientations of the particles (by the Brownian molecular motion) as well as all possible distances of the scattering centers

$$P(\vartheta) = \frac{1}{N} \sum_{i=1}^N \sum_{j=1}^N \frac{\sin(qr_{ij})}{qr_{ij}}$$

N = Total number of scattering centers

r = Distance of the scattering centers



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The fundamental equation of light scattering is corrected by the above mentioned **Debye interference theory** to

$$\frac{Kc}{R_\theta} = \frac{1}{MP(\theta)}$$

The scattering actuator $P(q)$ varies with the radius of gyration R_g according to

$$P(\theta) = 1 - (4\pi/\lambda_0 \sin \theta/2)^2 \frac{\langle R^2 \rangle}{3} + \dots$$


This results in the fundamental equation to

$$\frac{Kc}{R_\theta} = \frac{1}{M} \left(1 + (4\pi/\lambda_0 \sin \theta/2)^2 \frac{\langle R^2 \rangle}{3} \right)$$

From the plot of Kc/R versus q respectively q^2 with extrapolation to $q = 0$, the gyration radius R_g is obtained from the slope of the extrapolation straight line, derived from the extrapolations of the concentrations to $c \rightarrow 0$



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System





Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

Thus, the characteristic course of the angular dependent scattering intensity provides information about the shape and size of the scattering particle. Examples of form factor models are (among others):

Monodisperse Coil

$$P_z(u = qR_g) = \left(\frac{2}{u^4} (\exp(-u^2) - 1 + u^2) \right)^2$$

Polydisperse Coil

$$P_z(u = qR_g) = \left(1 + \frac{u^2}{3} \right)^{-1}$$

Thin Rod

$$P_z(u = qL = 12qR_g) = \frac{2}{u} \int_0^{\frac{u}{2}} \frac{\sin(u)}{u} du - \left(\frac{2}{u} \sin\left(\frac{u}{2}\right) \right)^2$$

L length of the rod,

Hard Sphere

$$P_z(u = qR_s = \sqrt{\frac{5}{3}} qR_g) = \left(\frac{3}{u^3} (\sin(u) - u \cos(u)) \right)^2$$

R_s radius of the sphere,



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

2.0. Performing a Static Light Scattering measurement

Sample Preparation

- preferably monodisperse particles
- determine an appropriate solvent (organic / aqueous, etc.)
- an accurate sample weight and volume determination for exact sample concentration
- determine/measure the refractive index of the solvent
- determine/measure the refractive index increment (dn/dc) of the particle in the corresponding solvent -> measurement with a differential refractometer
- perfect solution (of particles) required -> no loss of concentration
- filtration of the sample (PTFE filter for organic solvents / CA filter for aqueous solvents), use filter with appropriate pore size
- Preparation of concentration series required (at least for determination of molecular weight)

|



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System


2.1. Performing a SLS measurement using the ALV-Correlator Software



1.) Select the angular dependent measurement mode via **Quick Set**



2.) Click on **Dynamic/Static Light Scattering** and subsequently select **General DLS and SLS**





Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

3.) Angular Dependent DLS&SLS measurement mode

Angular Dependent Dynamic and Static Light Scattering

Angular Dependent Dynamic and Static Light Scattering

STEP 1 Measurement Control

Min. Angle: Max. Angle: Angular Step:

Runs: Duration[s]: **Accept this Set →**

Remeasure if dRate exceeds %

Sample **Setup** **Cumulant** **Atten.**

	Lower	Upper	Step	Runs
Range 1				
Range 2				
Range 3				
Range 4				
Range 5				

STEP 2 Static & Dynamic LS Options

Measure and correct for detector dark counts

Use Toluene Rayleigh Ratios as standard
... presuming a Refractive Index of

Compute Rayleigh Ratio for VV-Scattering
Instead use this RR : [1/cm]

Correct for index matching fluids refractive index
Refractive Index =

Compute and store Cumulants Analysis
Data cutoff at : % of maximum

Generate measurement protocol

STEP 3 Measurement Type

Standard Solvent Solution

Standard File : 

Solvent File : 

Solution File : 

Measure without solvent information/file

Append solution data to solution file

Create Angular Dependent Report

Sample Concentration : [g/ml]


Sample dn/dc : [ml/g]

Start Measurement **Cancel** **Help**



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

4.) Enter the favored measurement sequence in the **Measurement Control** section



	Lower	Upper	Step	Runs
Range 1				
Range 2				
Range 3				
Range 4				
Range 5				

- A) enter the starting angle (**Min. Angle**) and the final angle (**Max. Angle**) as well as the angular steps (in °)
- B) enter the number of runs (**Runs**) at each measurement angle, as well as the duration of each run (**Duration**)
- C) transfer the settings with **Accept this Set** to the measurement sequence chart
- D) if required use a „dust-filter“ function (standard deviation of the mean intensities) via **Remeasure if dRate exceeds ... %**
- E) enter the sample parameters (refractive index and viscosity of the solvent) or select an existing solvent from the schedule of the **Sample** menu
- F) select the correlator settings (Auto/Cross Correlation, etc.) in the **Setup** menu



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

5.) Setting the measurement parameters in section *Static & Dynamic LS Options*

STEP 2 Static & Dynamic LS Options

Measure and correct for detector dark counts **G**

Use Toluene Rayleigh Ratios as standard **H**

... presuming a Refractive Index of

Compute Rayleigh Ratio for VV-Scattering

Instead use this RR : [1/cm]

Correct for index matching fluids refractive index

Refractive Index =

Compute and store Cumulants Analysis **I**

Data cutoff at : % of maximum


Generate measurement protocol **K**

- G) tick the **Measure and correct for detector dark counts** box for the prior measurement of the detector's dark counts
- H) select the standard value for the **Rayleigh Ratios** of toluene
- I) activate the **Compute and store Cumulants-Analysis** function in order to determine the hydrodynamic radius R_H derived from the Cumulants Fit of the simultaneous DLS measurement
- K) select **Generate measurement protocol** to get the measurement data/graph indicated in the *Ratio/Scattering Angle Display*



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

6.) Selecting the measurement mode using the **Measurement Type** menu



L) Select the **Standard** mode for the standard measurement of toluene and save the data as Standard.TOL file (or use a different file name) respectively use a file of an already existing toluene measurement as Standard File



Click on **Start Measurement** to start the measurement !



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

STEP 3 Measurement Type M

Standard
 Solvent
 Solution

Standard File : STANDARD.TOL 

Solvent File : SOLVENT.TOL 

Solution File : 

Measure without solvent information/file
 Append solution data to solution file
 Create Angular Dependent Report

Sample Concentration : [g/ml]

Sample dn/dc : [ml/g]


- M) Select the **Solvent** mode for the measurement of the (pure) solvent used for the preparation of the sample solution and save the data as *Solvent.TOL* file respectively use a file of an already existing solvent measurement as Solvent File

Click on **Start Measurement** for starting the measurement !

Please note: The prior measurement of the (pure) solvent serves to determine the scattering intensity contribution of the solvent, in order to subtract it later on from the total intensity to get the excess scattering of the sample. The measurement of the solvent's intensity is only useful in those cases, where the scattering intensity of the solvent has already a significant amount of the total intensity. In any other cases this step can be skipped.



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



STEP 3 Measurement Type

Standard
 Solvent
 Solution N

Standard File : STANDARD.TOL

Solvent File : SOLVENT.TOL

Solution File : Test.STA

Measure without solvent information/file O
 Append solution data to solution file P, Q
 Create Angular Dependent Report

Sample Concentration : 1e-3 [g/ml]

Sample dn/dc : 0.1 [ml/g]

- N) for the measurement of the sample solution select the measurement mode **Solution** and save the corresponding data as *.STA data file (Solution file)
- O) In case no *Solvent File* was priorly measured (or loaded) tick the **Measure without solvent information/file** box – the scattering intensity of the solvent will be assumed as zero.
- P) Select **Append solution data to solution file** in case several sample concentrations will be measured and should be finally safed in one measurement file
- Q) **Create Angular Dependent Report** creates a report of angular dependency including distribution functions



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

STEP 3 Measurement Type

Standard
 Solvent
 Solution

Standard File : STANDARD.TOL 

Solvent File : SOLVENT.TOL 

Solution File : Test.STA 

Measure without solvent information/file
 Append solution data to solution file
 Create Angular Dependent Report

Sample Concentration : [g/ml] **R**

Sample dn/dc : [ml/g] **S**

- R) enter the value of the sample concentration (g/ml) in the **Sample Concentration** section
- S) enter the value of the refractive index increment (ml/g) of the sample particles (for the corresponding solvent used) in the **Sample dn/dc** section

Click on **Start Measurement** to start the measurement !




Statische Lichtstreuung mit dem ALV/CGS-3 Kompakt Goniometer System

Values of the refractive index increment $\frac{dn}{dC}$ for typical sample/particle classes dissolved in corresponding solvents


Sample/Particle class	Solvent	dn/dc [mL/g]
Lysozyme	Water /water based buffer	0,185
Phospholipide (Liposome)	Water	0,16
Polystyrene	cis-Decaline	0,12
Polystyrene	Toluene	0,08 – 0.11
Polystyrene	Tetrahydrofurane (THF)	0,18 – 0,19
Polymethylemethacrylate (PMMA)	Tetrahydrofurane (THF)	0,09
Polymethylemethacrylate (PMMA)	Toluene	0.01-0.02
Polyvinylchloride (PVC)	Tetrahydrofurane (THF)	0,10
Polyethyleneglycol	Water	0,13



3.0. Performing a SLS measurement using the ALV/Fit & Plot Software



Open the *.STA file of the measured sample via **File -> Load New File**



The screenshot shows the ALV software interface. The main window is titled "ALV / Static & Dynamic FIT and PLOT" and has a menu bar with "File", "Edit", "Window", and "Help". The "File" menu is circled in red. Below the menu bar are icons for file operations. A dialog box titled "NoName" is open, displaying "Loading File ! Please Wait !" and an "Abort Loading" button. To the right, a file selection dialog titled "Öffnen" is open, showing a list of files in the "ALVStat" folder:

Name	Änderungsdatum	Ty
Ps_std1.sta	31.10.2008 19:50	S1
Ps_std2.sta	24.03.1998 13:02	S1

The "Dateiname:" field is empty, and the "Dateityp:" is set to "ALV5000 [stat.-dyn.] (*.sta;*.dat)".



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



ALV-Stat: ALV-5000 Data

Fixed Parameter	
No. of Data Record	2
Name of standard sample	"Samplename needed."
Name of solvent sample	"Samplename needed."
Temperature/K	295.904

Adjustable Parameters	
n-Standard	0.00001
n_{Solvent}	1.49600
$\frac{dn}{dc} / 1/(g/dm^3)$	1.06000E-04
Concentration / g/dm^3	6.70000
Wavelength / nm	532.00
Rayleigh Ratio RR / 1/m	2.01700E-03
Viscosity / mPas (cP)	5.66330E-01
Sample name	ps sample1

dm³ = litre

Continue without any change

Abort

? Hilfe

(1) Prompt again after reading the next concentration

(2) Perform for all records/concentration do NOT prompt again. This is relevant only in case of concentration value itself not has been edited here. Otherwise program will turn to option

Change

Attention ! Make a copy of the original data file before pressing "Change", because recalculated data immediately overwrite the original data in the file !



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The *File Parameter Menu* shows a summary of the fixed respectively modifiable measurement parameters of the sample measurement – fixed (unchangeable) parameters are:

- selected *Standard.TOL* file
- selected *Solvent.TOL* file


Modifiable parameters (in retrospect) are:

- refractive index *Standard* (typically toluene)
- refractive index solvent (for DLS)
- refractive index increment dn/dc of sample particles for corresponding solvent
- sample concentration
- wavelength of used light source
- Rayleigh Ratio of the standard (typically toluene)
- viscosity of solvent (for DLS)
- sample name

The sample measurement file can be opened either without any change or previous used parameters will be changed (1) individually for each concentration measurement or (2) for all subsequent concentration measurements (of the same sample)



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System




The file data is indicated in the preset Zimm-Plot (Kc/R versus q^2+kc), by using the spread factor k the q -functions of each measured concentrations can be stretched in order to get the functions more properly separated.



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The appropriate Plot-Typen can be selected via the **Plot-Type** function; there are the following Plot-types for the evaluation of the Static Light Scattering measurements available:



Zimm-Plot (Kc/R vs. q^2+kc)

- Suitable for mono-disperse particles with particle size in the range of 20-50nm

The Zimm-Plot is a graphical evaluation method for scattering functions under consideration of the polydispersity of the macromolecules. Hereby the weight averaged molecular weight M_w can be derived from the scattering function Kc/R by double extrapolation, using the extrapolations $c \rightarrow 0$ (condition without any particle interactions) and $\theta \rightarrow 0$ (condition without scattering interference System). To get these functions the intensity I of the scattered light for different sample concentrations (same sample material) will be measured at different scattering angle positions (respectively q -vectors) as an angular dependent measurement. The measured intensities are referred to the priorly measured *Standard* (typically toluene) and actually expressed as Reduced Scattering Intensity I_{red} .

The obtained intensity values will be converted into ratio values and subsequently transferred in a plot, whose ordinate ($K \cdot c/R$) and abscissa ($q^2 + k \cdot c$) are depicted. After that the data points, which were measured at the same concentration respectively at the same scattering angle, will be interconnected.



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

The file data is indicated in the preset Zimm-Plot (Kc/R versus q^2+kc), by using the spread factor k the q -functions of each measured concentrations can be stretched in order to get the functions more properly separated.

The lines joining these points are ideally straight lines (first order analysis). The Zimm-Extrapolation will be derived from extending the obtained function lines towards the axis and subsequently determining the appropriately extrapolated data point of the corresponding abscissa value of q^2 respectively $K \cdot c$.

The hereby extrapolated values will be connected by a corresponding straight line and extended to the intersection point of the ordinate. The intercept value at the ordinate (y-axis) gives the inverse weight-average molecular weight ($1/M_w$). Additionally the Zimm-diagram provides the Radius of Gyration, derived from the slope of the extrapolation line $c \rightarrow 0$ respectively the second Virialcoefficient from the slope of the extrapolation line $q \rightarrow 0$.

Berry-Plot ($\sqrt{Kc/R}$ vs. q^2+kc)


- appropriate for samples with higher molecular weight , typ. diameter > 50nm

Guinier-Plot ($\ln Kc/R$ vs. q^2+kc)

- appropriate for small angles and large particles, $qR_g < 1$



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



$M_w(c): 9.286e+05 \text{ g/mol}$


$M_w(q^2): 9.286e+05 \text{ g/mol}$

$A_2: 5.522e-07 \text{ mol dm}^3/g^2$

$R_g: 4.357e+01 \text{ nm}$



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



via **Calculation** -> **Results** the result spread sheet can be obtained, in which molecular weight M_w , static structure factor $\langle S^2 \rangle$ and radius of gyration R_g for each concentrations as well as the extrapolated ($c=0$) concentration are plotted.

ALV-Stat: Fit Results

4.51 10/10
 Static Light Scattering, results of Zimm-Plot, q^2 -dep.: 1, c -dep.: 1

File	Conc/(g/dm ³)	Mw(app)/(g/mol)	$\langle S^2 \rangle$ (app)/ μm^2	R_g
Conc.=0	0.0000	9.286e5	1.898e-3	
Ps_std1.sta	6.7000	1.167e5	1.819e-4	
Ps_std1.sta	4.6400	1.641e5	2.447e-4	
Ps_std1.sta	3.2400	2.155e5	4.333e-4	
Ps_std1.sta	1.9100	3.237e5	5.819e-4	
Ps_std1.sta	0.9700	4.393e5	8.691e-4	

Mw(c): 9.286e+05 g/mol ($\pm 10.2\%$) Mw(q^2): 9.286e+05 g/mol ($\pm 1.67\%$)
 A2: 5.522e-07 mol dm³/g² ($\pm 2.46\%$) Rg: 4.357e+01 nm ($\pm 2.05\%$)

Results to Clipboard OK



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System



Hard Sphere
with Radius R

$$R_g^2 = \frac{3}{5} R^2$$

Spherical Capsule
with outer radius R_o
and inner radius R_i

$$R_g^2 = \frac{3}{5} \frac{R_o^5 - R_i^5}{R_o^3 - R_i^3}$$

Cylinder
with radius R and
length h

$$R_g^2 = \frac{R^2}{2} + \frac{h^2}{12}$$

Hollow Cylinder
with outer radius R_o
and inner radius R_i
and length h

$$R_g^2 = \frac{R_o^2 + R_i^2}{2} + \frac{h^2}{12}$$

Ellipsoide
with semi-axis a, b and c

$$R_g^2 = \frac{a^2 + b^2 + c^2}{5}$$

Disc-shaped Particle
Flat disc with radius R

$$R_g^2 = \frac{R^2}{2}$$

**Dependency of Gyrationradius R
for typical particle geometrics**



Static Light Scattering using the ALV/CGS-3 Compact Goniometer System

Four different form factor models can be selected for the plot of the form factor vs. the scattering vector q :

